



RTD™ WOUND DRESSING INSTRUCTIONS FOR USE



DESCRIPTION:

RTD™ Wound Dressing is made of polyurethane/ polyether foam, Methylene Blue (0.25 mg/g), Gentian Violet (0.25 mg/g), and Silver Sodium Zirconium Phosphate (7 mg/g).

INDICATIONS:

The RTD™ Wound Dressing is designed for use with adults only and for the management of partial to full thickness wounds with light to heavy exudate, and has demonstrated effectiveness for pressure ulcers, venous stasis ulcers, and diabetic foot ulcers.

APPLICATION: Easy as 1-2-3

Measure and Cut the RTD™ Wound Dressing to the size of the wound. The RTD Wound Dressing can be rolled or layered and packed into a wound.

Apply the RTD™ Wound Dressing in direct contact with the wound bed. Direct application is necessary to obtain the best results.

Cover with a non-occlusive secondary dressing of choice. Leave in place for up to 3 days.

CONTRAINDICATIONS:

RTD™ Wound Dressings are not indicated for burns; pregnant or breastfeeding women; or children.

PRECAUTIONS:

In the event of any clinical infection, the RTD™ Wound Dressing is not intended to replace systemic therapy or the appropriate course of treatment for infections.

When choosing a wound dressing, clinicians should consider any potential adverse effects for any patient.

Stop using straight away if an allergic reaction is suspected.

Caution:

POM. Keep out of reach of children. This device may only be sold by or on the order of a registered physician.

RECOMMENDATIONS:

Wound Preparation: Before applying the RTD™ Wound Dressing, clean the wound with an appropriate cleaning agent in accordance with wound care guidelines .

Standard precautions should be taken to prevent the spread of infectious agents when dressing wounds.

Securing RTD™ Wound Dressing In Place:

Use gauze®, Kerlix®, Omnifix®, Mefix® or other non-occlusive adhesives or medical grade tape.

Guidance: For best results, the RTD™ Wound Dressing should be in direct contact with the wound bed and fit to the edges of the wound.

For moderate to heavy exuding wounds the 6.35mm thickness is recommended. For light exuding and/or superficial wounds the 3.18mm thickness is recommended.

Leave the RTD™ dressing in place for three days or remove if saturated. If the dressing becomes saturated before the desired dressing change, a second dressing may be placed on top for continued absorption. If The RTD™ dressing is difficult to remove, soak with saline or sterile water. A non-adherent layer is recommended to prevent sticking or in-growth at the later stages of wound healing.

For Use At Home: The RTD™ Wound Dressing should only be used under the direct supervision of a wound care professional.

STERILE R Rx Only



Do not use if package is opened or damaged.



See Instructions.



Single use only.



Keep in a dry place away from direct sunlight.



Avoid storing at temperatures higher than 44°C/ 112°F



Manufactured for Keneric Healthcare, LLC • 825 Watters Creek Blvd, Suite 250 • Allen, TX 75013, USA • +1 855-872-2013

RTDH-IFUP-20190226CE



Obelis s.a., Bld.
General Wahis 53,
1030 Brussels, Belgium

©2018 Keneric Healthcare