

RTD™

WOUND DRESSING

PATIENT CARE IS YOUR PRIORITY. IT'S OURS TOO.



A VERSATILE DRESSING FOR TREATING ACUTE AND CHRONIC WOUNDS

- Creates Optimal Environment for Wound Healing
- Light to Heavy Exudate
- Effective on Pressure Ulcers, Diabetic Foot Ulcers and Venous Stasis Ulcers
- Ready to Apply and Easy to Use
- Increase Clinician Efficiency
- Streamline Facility Expenditures




www.KenericHealthcare.com

FOLLOW US    

RTD™ WOUND DRESSING



RTD™ Wound Dressing is the universal quick-closing, moist, advanced oxygenation, capillary-suction-therapy (CST) and biofilm reducing dressing with bio-burden sequestration.



RTD™ is a proprietary highly absorbent antimicrobial foam with Methylene Blue (0.25mg/g), Gentian Violet (0.25mg/g) and Silver Sodium Zirconium Phosphate (7mg/g) tightly bound into the polymer matrix.



Designed for the management of partial to full thickness wounds with light to heavy exudate, RTD™ has demonstrated effectiveness for pressure ulcers, venous stasis ulcers, and diabetic foot ulcers.

RTD™ FEATURES AND BENEFITS

Interactive ingredients: Silver, Methylene Blue, Gentian Violet and CST foam

- Draws protein rich exudate away from wound
- Highly effective antimicrobial and antifungal properties
- Early and sustained antimicrobial protection
- Significant pain reduction reported
- Non-cytotoxic
- Cost effective
- One dressing for all your needs

RTD™ PROPRIETARY TECHNOLOGY

- Interactive ingredients that are tightly bound into the RTD™ PU foam
- Open cell structure provides powerful vertical wicking of exudate
- Absorbs 1,110 mg/cm²
- Holds 15x its weight in exudate
- Regulates wound bed moisture with little to no maceration

RTD™ PU FOAM IS SOFT, PLIABLE AND SURFACE CONFORMING

- Easy to Use - Ready to Apply - Cut to Size - No Need to Hydrate

0.64cm THICKNESS - LIGHT TO HEAVY EXUDATE

KHC#	SIZE	QTY	HCPCS
76020214	5cmx5cm	10/BOX	A6209
72040414	10cmx10cm	20/BOX	A6209
76040414	10cmx10cm	10/BOX	A6209

**CONTACT US FOR
COMPLIMENTARY PRODUCT
TRIAL SAMPLES.**

AEscobedo@KenericHC.com



KENERIC
HEALTHCARE

www.KenericHealthcare.com

FOLLOW US    